---

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The screenshot shows the RDE web interface. At the top left is the RDE logo. At the top right are links for 'DICE', 'ログアウト', and '利用方法'. Below the header is a blue bar with '仕掛データセット一覧' and a button for 'データ登録状況一覧'. The main content area is titled '検索条件' and contains the following elements:

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A red bracket on the left side of the search criteria section points to the input fields, indicating that they can be used for direct search.

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The screenshot shows a web interface for selecting a device. It is titled "装置選択" (Device Selection). Under "検索条件" (Search Conditions), there is a dropdown menu for "機関" (Institution) which is highlighted with a red box and a red circle containing the number 1. Below it are fields for "事業" (Project), "装置タイプ" (Device Type), and "ローカルID" (Local ID). At the bottom of the search section, there is a "検索" (Search) button highlighted with a red box and a red circle containing the number 1, and a "リセット" (Reset) button. Below the search section is the "検索結果" (Search Results) section, which contains a table with two rows. The first row is highlighted with a red box and a red circle containing the number 2. The table has columns for "装置(登録名)" (Device (Registered Name)), "メーカー名" (Manufacturer Name), "型別名・型番" (Model Name/Number), "機関" (Institution), "事業" (Project), "ローカルID" (Local ID), and "装置タイプ" (Device Type). The first row has a radio button selected, and the second row has a radio button unselected. Below the table, there is a "Showing 1 to 2 of 2 entries" message and navigation buttons: "First", "Previous", "Next", and "Last". The "Next" button is highlighted with a red box and a red circle containing the number 3.

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The screenshot shows the '送り状入力' (Shipping Label Input) screen for a 'TESTデータセット' (TEST Data Set). The form contains the following information:

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The screenshot shows a web interface for sample selection. At the top, a blue header bar contains the text "試料選択" (Sample Selection), highlighted with a red box and a pink circle containing the number 3. Below the header is a search bar with the placeholder text "試料名、タグ、説明、化学式・組成式・分子式などを入力してください。" (Please enter sample name, tag, description, chemical formula, composition formula, molecular formula, etc.) and a magnifying glass icon, highlighted with a red box and a pink circle containing the number 4. Below the search bar is a table with three columns: "試料名(ローカルID)" (Sample Name (Local ID)), "試料の説明" (Sample Description), and "化学式・組成式・分子式など" (Chemical formula, composition formula, molecular formula, etc.). The table contains several rows of data, with the first row highlighted by a red box and a pink circle containing the number 5. At the bottom of the table, there is a "選択" (Select) button highlighted with a red box and a pink circle containing the number 6. Below the table, there is a pagination bar showing "Showing 1 to 100 of 1188 entries" and navigation buttons for "First", "Previous", "1", "2", "3", "4", "5", "...", "12", "Next", and "Last".

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The screenshot shows a web interface for data registration. It features two main sections: '登録ファイル' (Registered File) and '添付ファイル' (Attachment File). Both sections have a dashed box for file upload and a 'ファイル選択' (File Selection) button. A red box highlights the '登録ファイル' section, and a red circle with the number '1' is placed over the 'ファイル選択' button. At the bottom, there are '内容確認' (Check Content) and 'キャンセル' (Cancel) buttons.

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The screenshot shows a web interface for data registration. It features two main sections: '登録ファイル' (Registered Files) and '添付ファイル' (Attachments). Each section contains a dashed box for file upload and a 'ファイル選択' (File Selection) button. At the bottom, there are two buttons: '内容確認' (Content Confirmation) and 'キャンセル' (Cancel). The '内容確認' button is highlighted with a red box, and a red circle with the number '1' is placed next to it, indicating the step to click this button after data entry is complete.

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The screenshot shows a web form for data registration. At the bottom of the form, there are four buttons: '登録開始' (Register Start), '連続登録' (Continuous Registration), '修正' (Correction), and 'キャンセル' (Cancel). The '登録開始' button is highlighted with a red rectangular box. A red arrow points from this box to a larger, magnified view of the '実験データ' (Experiment Data) section on the right.

This is a magnified view of the '実験データ' (Experiment Data) section. It contains a table with the following data:

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Below the table, there are four buttons: '登録開始' (Register Start), '連続登録' (Continuous Registration), '修正' (Correction), and 'キャンセル' (Cancel). The '登録開始' button is highlighted with a red rectangular box. A red arrow points from this box to the '登録開始' button in the main form screenshot on the left.

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The screenshot displays a web interface for data management. At the top, there is a blue header with the text 'データ一覧:' (Data List). Below the header, there are two buttons: 'データセット詳細' (Data Set Details) and 'データセット一覧へ戻る' (Return to Data Set List). The main content area shows a list of datasets. The first entry is 'sample\_data.TXT', which is highlighted with a red rectangular box. Below the dataset name, there is a preview image showing the word 'SAMPLE' on a document. To the right of the dataset name, the number '1' is displayed. Below the preview, there is a section for 'データ所有者(所属)' (Data Owner) and '資料名(ローカルID)' (Material Name (Local ID)), with the value '操作体験会用サンプル' (Sample for Operation Experience). Other details include '説明' (Description) '操作体験会' (Operation Experience), 'タグ/シノミー' (Tag/Synonym) 'ワクシノミー' (Vaccinomy), and '登録日時' (Registration Date) '2022-12-05 12:09:26 JST'. At the bottom of the list, it says 'Showing 1 to 1 of 1 entries'. Navigation buttons for 'First', 'Previous', 'Next', and 'Last' are visible at the bottom right.

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The screenshot shows a web interface for data details. At the top, there is a blue header with the text "データ詳細:". Below the header, there are three buttons: "データ削除" (Data Delete), "データダウンロード" (Data Download), and "データ一覧へ戻る" (Return to Data List). A red box highlights three tabs: "概要" (Overview), "ファイル" (Files), and "添付ファイル" (Attachments). Below the tabs, there is a "送り状編集" (Edit Shipping Label) button. The main content area displays a graph with "Absorbance (a.u.)" on the y-axis and "Wavelength (nm)" on the x-axis. The graph shows a single sharp peak. A large "SAMPLE" watermark is overlaid on the graph. Below the graph, the filename "sample\_data\_extract.png" is visible. At the bottom, there is a section for "サムネイルをクリックすると画像が表示されます。" (Click the thumbnail to display the image.) with two thumbnails labeled "sample\_data\_ext" and "fact.png".

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The screenshot shows a web interface for 'Data Details'. At the top, there are three buttons: 'データ削除' (Delete Data), 'データダウンロード' (Download Data), and 'データ一覧へ戻る' (Return to Data List). Below these is a navigation bar with three tabs: '概要' (Summary), 'ファイル' (Files), and '添付ファイル' (Attachments). The '概要' tab is highlighted with a red box. To the right of the tabs is a '送り状編集' (Edit Invoice) button. The main content area features a large graph with 'Absorbance (a.u.)' on the y-axis and 'Wavelength (nm)' on the x-axis. A prominent 'SAMPLE' watermark is overlaid on the graph. Below the graph, there is a note: 'サムネイルをクリックすると画像が表示されます。' (Click the thumbnail to display the image.) and two thumbnails labeled 'sample\_data\_ext\_rect.png'. At the bottom of the page, there is a section labeled 'メタ情報' (Meta-information).

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The screenshot displays a web interface for data management. At the top, there is a blue header with the text 'データ一覧:' and two buttons: 'データセット詳細' and 'データセット一覧へ戻る'. Below the header, there is a search bar with the text 'ギャラリー表示' and a search button '> 検索条件'. On the right side, there are filters for '登録日時' (dropdown), '状態' (radio buttons), and '種類' (radio buttons). The main content area shows a list of data entries. The first entry is 'sample\_data.TXT' with a small icon to its right. This entry is highlighted with a red rectangular box. Below the entry name is a thumbnail image with the word 'SAMPLE' overlaid. Underneath the thumbnail, there is a section for 'データ所有者(所属)' followed by 'ID' and 'ID' labels, and a '登録日時' field showing '2022-12-05 12:00:29 JST'. At the bottom of the page, there is a pagination control showing 'Showing 1 to 1 of 1 entries' and buttons for 'First', 'Previous', 'Next', and 'Last'.

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The screenshot displays the '送り状編集' (Shipping Label Edit) interface. At the top, a blue header bar contains the text '送り状編集'. Below this, a white panel shows metadata: 'データセットタイプ' (Data Set Type) as '構造解析・リファレンス型' (Structural Analysis/Reference Type), '更新日時' (Update Date) as '2023-02-08 16:13:41 JST', and '対象装置' (Target Device) as '課題番号/課題名' (Task Number/Task Name). A '基本情報' (Basic Information) section follows, with '記入年月日' (Entry Date) as '2023-02-08 JST' and '装置' (Device) as '装置' (Device). The bottom section, '添付ファイル' (Attachments), features a dashed box with the instruction '添付ファイルをドラッグアンドドロップしてください。' (Please drag and drop the attachment file here.) and a blue 'ファイル選択' (File Selection) button. At the very bottom, there are '保存' (Save) and 'キャンセル' (Cancel) buttons.

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The screenshot displays a web interface for data management. At the top, there is a blue header with the text 'データ一覧:' (Data List). Below the header, there are two buttons: 'データセット詳細' (Data Set Details) and 'データセット一覧へ戻る' (Return to Data Set List). A search bar is present with the text '検索条件' (Search Condition). On the right side, there are filters for '登録日時' (Registration Date), '公開' (Public), and '閲覧' (View). The main content area shows a list of datasets. The first entry is 'sample\_data.TXT', which is highlighted with a red rectangular box. Below the dataset name, there is a preview image showing a graph with the word 'SAMPLE' overlaid. Further down, there are sections for 'データ所有者(所属)' (Data Owner), '試料名(ローカルID)' (Sample Name), '説明' (Description), and 'タグ/ノミー' (Tag/Name). At the bottom, there is a pagination control showing 'Showing 1 to 1 of 1 entries' and buttons for 'First', 'Previous', 'Next', and 'Last'.

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The image shows a web interface for sample registration. On the left, a form titled '試料情報' (Sample Information) contains several input fields. The '試料名(コールID)' (Sample Name (Call ID)) field is highlighted with a red box and contains the text 'test-1'. A red arrow points from this field to a larger, detailed view of the same field on the right. In this detailed view, the text 'test-1' has been replaced with 'new sample 2023'. The '新規登録' (New Registration) checkbox is checked. At the bottom of the main form, the '保存' (Save) button is highlighted with a red box.

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